152321.0200

Stereomikroskop BVF

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Stereo Microscope BVF

Operating instructions 5

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Mode d'emploi 8

Microscopio stereoscopico BVF

Manuale d'uso 11

Microscopio estéreo BVF

Manual de instrucciones 14

BVF sztereomikroszkóp

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Mikroskop stereoskopowy BVF

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Stereomicroscop BVF

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Stereomikroskop BVF

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152321.0200

STEREOMIKROSKOP BVF

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Einleitung

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Operating instructions

Introduction

In purchasing a NERIOX stereo microscope, you have chosen a quality product. It is suitable for use in schools and laboratories as well as in industry.

The Stereo Microscope consists of two separate microscope tubes, which are combined into one unit to allow simultaneous focusing on the object. Each tube has a prism, an achromatic lens and a wide field eyepiece. This creates a large, flat field of view. Both eyes view the object at a slightly different angle, which creates a stereoscopic effect.

With normal use, the maintenance requirement is reduced to a minimum.

These operating instructions provide information on the design, handling and maintenance of the microscope.

Design of the Stereo Microscope

Below you will find the designation of the individual components, which are also shown in the illustration:

1. Eyecup
2. Dioptre adjustment
3. Lens revolver
4. Switch for lighting
5. Object clamps
6. Object plate
7. Stand base with lighting
8. Incident illumination
9. Focusing (2 rotary knobs)
10. Prism housing
11. Wide field eyepiece
12. Rotary knob light

Operating principle of the individual components

The Stereo Microscope consists of a stand base, a vertical stand column and a head with focusing device. Always hold the microscope by the stand when you want to move it.

Tubes

The two tubes have a tilt angle of 45° and can be rotated by 360° if necessary.

Optical specifications

The Stereo Microscope has a pair of 10x wide field eyepieces and two pairs (non-interchangeable) of achromatic 2x and 4x lenses in a rotary lens revolver. The total magnification of the microscope is calculated by multiplying the magnification of the eyepieces by the magnification of the lenses. The eyepieces are each fixed in the tube with a small screw.

Object stage

The object stage has integrated transmitted light illumination, a transparent object plate and two object clamps. A black and white object plate is also provided.

Attention:

When using the black and white object plate, the transmitted light illumination must be switched off!

Focusing

Focusing is used to bring the image of the object into focus. The working distance (distance between the underside of the lens and the top of the object) is approx. 75mm.

Lens revolver

Two pairs of achromatic lenses are integrated in the rotary lens revolver. A different magnification is set by rotating the lens revolver by 180°. Turn the lens revolver until you hear a click.

Illumination

This model is equipped with LED incident light and transmitted light illumination

- LED: Incident light illumination
- LED: Transmitted light illumination

Both lights can be used simultaneously (switch for lighting).

Rechargeable batteries are optionally available for the LED lighting on these models. The charging time is approximately 15 hours. The capacity of the charged batteries is sufficient for about 35 hours of operation.

Batteries: 3 nickel-metal hydride batteries, type AA, 1.2 volt 2000 mA. → 155060.0220

Preparing the Stereo Microscope for use

Remove the microscope from its packaging and place it on a stable table. Connect the power supply and insert the mains plug into a power socket. Switch on the lighting with the switch for lighting. Then slide the eyecups over the eyepieces. Now sit comfortably in front of the microscope and look through the eyepieces.

Working with the Stereo Microscope

Place an object on the object plate and set the ideal height for the object with the height adjustment mechanism on the stand column.

Then use the focusing mechanism to focus the image.

- Adjust the pupil distance by moving the prism housings towards or away from each other until you see a round image with both eyes.

- Normally, both eyepiece tubes should have the same length. To adjust for any difference between the two eyes, please proceed as follows:

Take a flat object and look into the right eyepiece with your right eye only. Now use the focusing mechanism to focus the image.

Now look through the left eyepiece with your left eye only and use the dioptre adjustment mechanism on the left tube to get a sharp image.

- The models with illumination also have a black and white object plate that is easy to replace. The black side can be used to observe objects that are better displayed against a dark background.

Maintenance and cleaning

Always cover the microscope with a dust cover after use. Keep the eyepieces mounted in the microscope tube at all times, to prevent dust from getting into the device.

Cleaning the optics

If the lenses of the eyepiece are dirty, they can be wiped with a piece of lens paper (circular movements). If this does not bring an improvement, moisten the lens paper with a drop of xylene or alcohol. Never apply xylene or alcohol directly on the lenses!

If dirt is clearly visible in the field of view, this is located on the lowest lens of the eyepiece. The eyepieces can be removed from the tube with the aid of a Phillips screwdriver (size 00). Now clean the outside of the lens. If dust is still visible, make certain that the dust is really in the eyepiece by rotating the eyepiece. If this is definitely the case, carefully unscrew the bottom lens from the housing and clean it.

Warning:

Cloths with plastic fibres can damage the coating of the lenses!

Cleaning the stand and toothed rack

Dust can be removed with a brush. If necessary, the toothed rack and the guide of the focusing mechanism can be treated with a little grease or oil.

Procedure for illumination problems (incident light / transmitted light)

The built-in LEDs have a very long service life. If there is a problem with the illumination, please send the microscope to our Service Center → Servicecenter@brw.ch.

Warning:

LEDs or electronic parts may only be repaired or replaced by qualified personnel! In the event of non-compliance, any liability will be refused.

Charging the batteries for battery operation

Before using the device for the first time in the battery mode, the batteries must be charged. Connect the charger cable to the connector on the back of the stand base and connect the charger to a power socket.

Replacing the batteries

Warning:

Always pull out the mains plug before replacing the batteries!

- Remove the cover of the battery compartment from the base plate of the stand.
- Remove the old batteries, insert the new batteries (type AA) in the correct direction and put the cover back on.

Repairs

If your microscope needs to be repaired, please contact our Service Center → Servicecenter@brw.ch

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